


Oxygen Concentrator Unit, MST(M2)-109U


Specification

Model Name	MST(M2)-109U
Oxygen Flow	10LPM
Oxygen Purity	Above 90%
Oxygen Pressure	0.95~1.0bar
Voltage(Frequency)	230V/50Hz
Power Consumption	Below 650W
Dimension(W x D x H)	230 x 490 x 545mm
Weight	Approx.20kg



Oxygen Tank + Booster Compressor

- Oxygen can be supplied individually to each site that needs oxygen supply
- The unit can be used by connecting the oxygen tank and the booster compressor




RAK-U06M2E

- RAK-U06M2E is Oxus's Product
- Units can be individually turned on/off depending on the amount of oxygen required.



Oxygen Concentrator Unit, MST(M2)-109U

- Oxygen Unit
- Each unit produces oxygen



Oxygen Tank + Booster Compressor

- Oxygen can be supplied individually to each site that needs oxygen supply
- The unit can be used by connecting the oxygen tank and the booster compressor



RAK-U06M2E

- RAK-U06M2E is Oxus's Product
- Units can be individually turned on/off depending on the amount of oxygen required.

Specification

Model Name	Oxygen Flow	Oxygen Purity	Oxygen Pressure	Voltage(Frequency)	Power Consumption	Dimension(W x D x H)	Weight
MST(M2)-109U	10LPM	Above 90%	0.95~1.0bar	230V/50Hz	below 650W	230 x 490 x 545mm	App.15kg